1010670

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Version 1.0

ATTENTION!

- Not suitable for children under 3 years
- Please read the instructions for use before the first use
- Use only under supervision of adults
- Use only on non-slip level surface
- Use GymTop MODO barefoot only or with non-slip, closed shoes

General information

- Make sure that sufficient practice room is available and no edges, corners or objects represent a risk of injury.
- During your first exercises on the GymTop MODO you probably would like to hold on or support. Please pay attention to suitable, reliable device aids or ask for assistance in partnership.
- GymTop MODO –like the classic therapy top– is the ideal training device, to train and improve versatilely general movement safety, coordination, response capacity, body posture, body perception and balance.
- Do not use GymTop MODO after intensive physical exertion, in case of pain, tiredness and physical or mental indisposition.
- Requirements for training with open eyes:
Stand on a level surface one after the other for about 15 sec. each on your right and left leg, without touching the floor with the other.
- Requirements for training with closed eyes:
Now perform this exercise with your eyes closed.

- Weight load up to 120 kg
- Material: 90% recyclable PP, 10% organic pigments – sturdy, abrasion-proof, free of harmful substances, washable, for inside and outside
- Dimensions : 40 cm diameter 8,5 cm height
- Weight: 364 g

Exercises

- Stand upright with both feet and slightly bent knees on the device und balance yourself GymTop MODO.
- Swing on the device in alternation forward and backward, to the left and to the right, whereby the upper body should remain upright. The device should be controlled by leg and hip movements.
- Stand upright with both feet and slightly bent knees on the device and try, as long as possible to keep the balance and remain in this starting position.
- Try to maintain balance when standing on one leg, practice in rotation with your right and left leg.
- In one-legged stand stretch on leg each in alternation laterally, to the front and to the back.
- If you feel relatively safe, you may try out some of these exercises with closed eyes. If necessary first with assistance.
- Make slight knee bends from a standing position or crouch down carefully. Thereby keep the upper body upright.
- Clap your hands while standing in front and behind the body together in (rapid) alternation.
- Get out of the two-legged stand on tiptoe. Try this exercise in the binding stand.
- Throw a ball back and forth for example with a partner or against a wall and catch it.

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©Jakobs GmbH